

OptiSmooth Mitro Catheter with 4 drainage eyes for Single Use Only

LATEX FREE

Thank you for choosing the OptiSmooth Mitro catheter with 4 drainage eyes. It has been specially developed from many years of listening to people who catheterise. Consulting experienced healthcare professionals we think it is one of the most advanced intermittent catheters available today.

This product is made of medical grade, high quality thermosensitive PVC and is intended for Single Use Only.

The OptiSmooth Mitro catheter incorporates 4 staggered drainage eyes developed to improve urine drainage.

It is available in two lengths a 40cm and 50cm.

Product type or specification

40cm	OSM4010	OSM4012	OSM4014	OSM4016	OSM4018
50cm	OSM5010	OSM5012	OSM5014	OSM5016	OSM5018

Intended Use

This product is to be used for intermittent self-catherisation (ISC) purposes and urine drainage purposes only

Precautions

- This product should be used solely for intermittent catheterisation , upon the recommendation of a qualified healthcare professional
- For single use only. Do not reuse.
- Do not use if the package or product has accidentally been damaged or contaminated.
- Use immediately after opening. After using please dispose of the catheter by placing the used catheter back into the opened packaging and dispose the pack safely and responsibly

Recommended method of use

1. Begin by washing your hands
2. Open the packet and take out the catheter using the funnel. Coat the catheter with a water soluble gel lubrication.
3. Position the OptiSmooth Mitro catheter at the urethral or Stoma (Opening/channel) with one hand; using the other hand to insert the lubricated catheter gently down into the urethra or passage until urine starts to flow.

You can drain the urine straight into the toilet or a receptacle or you can connect the OptiSmooth funnel to a urine drainage bag if preferred.

4. When the bladder is completely empty, slowly remove the catheter and dispose of it by placing it in the empty product pack. Please dispose of the product in an appropriate waste bin

Storage recommendation

Do not store at extreme temperatures and avoid humidity.

Avoid direct sunlight

OptiSmooth[®] Mitro

Optimum smoothness Optimum comfort

LATEX FREE

OptiSmooth[®] Mitro

Optimum smoothness Optimum comfort

OptiSmooth Mitro Catheter with 4 drainage eyes
for Single Use Only

 **Zhejiang Sungood Technology Co.,Ltd**
Add: No.39, 3rd Road Baiyang Street Economic Development
Zone Hangzhou, Zhejiang, China
Tel: 0086-571-28969218 Fax: 0086-571-28969226
http://www.sungoodmed.com.cn
E-mail: *****@sungoodmed.com.cn

 **Curas Ltd.**
289 Conston Road, Bristol, BS3 1RT, England
Tel: +44 7796 338585
E-mail: info@curas.com

CE 0197

STERILE EO



Vesica Urology Ltd

OptiSmooth[®] Mitro

Optimum smoothness Optimum comfort

INSTRUCTIONS FOR USE

Vesica Urology Ltd

PO Box 4145
Slough
Berkshire
SL1 0QA

Phone: +44 (0) 1753 445087

Fax: +44 (0) 1753 447446

enquiries@vesicaurology.com

www.vesicaurology.com