



OptiSmooth®

Optimum smoothness Optimum comfort

Hydrophilic Coated Intermittent Catheter
Requires additional water to be added
(For single use only)









Vesica Urology Ltd

OptiSmooth Key Features






- Supremely smooth and stable advanced hydrophilic coating which activates within 30 seconds contact with water. Unlike many other hydrophilic catheters the coating stays on the catheter and also resists drying out providing optimum comfort on insertion and removal.
- A new short Female Mini 12cm catheter is now available providing optimum discretion.
- Two finger holes to aid opening of the packaging. Just peel apart in opposite direction to gain access to the catheter.
- Flexible funnel connects securely to an external collection device such as a night drainage bag. The funnel has also been developed to secure firmly to a syringe for instillation procedures.
- Sticky Pad allows the OptiSmooth catheter to be secured to a clean surface allowing you time to prepare.
- Foil inner packaging protects the OptiSmooth catheters to ensure they reach you in optimum condition.
- Meaningful product codes OS for OptiSmooth the next two numbers represent the catheter length and the final two numbers are the carrier size (catheter thickness).

OptiSmooth Product Range






Male / Mitrofanoff (40cm length)

| Size | Funnel Colour | Product Code | PIP Code |
|------|--|--------------|----------|
| CH08 |  | OS4008 | 354-9573 |
| CH10 |  | OS4010 | 354-9581 |
| CH12 |  | OS4012 | 354-9599 |
| CH14 |  | OS4014 | 354-9607 |
| CH16 |  | OS4016 | 354-9615 |
| CH18 |  | OS4018 | 354-9623 |




Female (18cm length)


| Size | Funnel Colour | Product Code | PIP Code |
|------|--|--------------|----------|
| CH08 |  | OS1808 | 354-9631 |
| CH10 |  | OS1810 | 354-9649 |
| CH12 |  | OS1812 | 354-9656 |
| CH14 |  | OS1814 | 354-9664 |
| CH16 |  | OS1816 | 354-9672 |

Female Mini (12cm length)

| Size | Funnel Colour | Product Code | PIP Code |
|------|---|--------------|----------|
| CH08 |  | OS1208 | 354-9680 |
| CH10 |  | OS1210 | 354-9698 |
| CH12 |  | OS1212 | 354-9714 |
| CH14 |  | OS1214 | 354-9706 |
| CH16 |  | OS1216 | 354-9722 |

Paediatric (30cm length)

| Size | Funnel Colour | Product Code | PIP Code |
|------|--|--------------|----------|
| CH06 |  | OS3006 | 354-9730 |
| CH08 |  | OS3008 | 354-9748 |
| CH10 |  | OS3010 | 354-9755 |



Vesica Urology Ltd
PO Box 4145
Slough
SL10QA

Phone : +44 (0) 1753 445087
Fax : +44 (0) 1753 447446
Mobile : +44 (0) 7948 172607
adrian@vesicaurology.com
www.vesicaurology.com

LATEX FREE

Vesica Urology Ltd

OptiSmooth® is a registered trademark
 OPTISMOOTH IS A UK UNREGISTERED DESIGN 2011