


## **OptiSmooth MD Meatal Dilator Key Features**

- Anatomically designed No Touch grip handle developed specially for a more controlled comfortable insertion and removal. The innovative grip has 3 unique indentations for two forefingers and a thumb to allow complete control of the dilator prior to insertion. Each indentation has raised micro areas to give further grip and control.
- Tapered tip helps with the insertion with meatal strictures female urethral tight spots and stenosed stomas.
- Supremely smooth and stable advanced hydrophilic coating which activates within 30 seconds of contact with water. Unlike many other hydrophilic dilators the coating stays on the dilator and also resists drying out, providing optimum comfort on insertion and removal. Requires water to be added to product pouch (tap water is fine).
- Two finger holes to aid opening of the packaging. Just peel apart in opposite direction to gain access to the dilator.
- Sticky Pad allows the OptiSmooth MD meatal dilator to be secured to a clean surface allowing you time to prepare.
- Meaningful product codes OSMDIL for OptiSmooth MD meatal dilator, the final two numbers are the charriere size (catheter thickness).
- FOR SINGLE USE ONLY

## **OptiSmooth MD Meatal Dilator Product Range**

(14cm length)

Size	Grip Colour	Product Code	PIP Code
CH08		OSMDIL08	409-6871
CH10		OSMDIL10	409-6889
CH12		OSMDIL12	409-6905
CH14		OSMDIL14	409-6913
CH16		OSMDIL16	409-6897
CH18		OSMDIL18	409-6988



Vesica Urology Ltd PO Box 4145 Slough SL10QA

Phone: +44 (0) 1753 445087 Fax: +44 (0) 1753 447446 Mobile: +44 (0) 7948 172607 adrian@vesicaurology.com www.vesicaurology.com

## **Vesica Ur**logy Ltd

OptiSmooth® is a registered trademark
OptiSmooth MD® is a registered trademark
OptiSmooth MD meatal dilator is a UK registered design
No's UK 6046935, UK 6046936 & UK 6046937
EUIPO Registered design No's 003820885-001, 003820885-0002
CN ZL2017304521638