

# OptiSmooth<sup>®</sup> Mitro<sup>®</sup> XL<sup>®</sup>

Optimum smoothness Optimum comfort

**Uncoated nelaton catheter with 4 unique oversized  
xtra large (XL) Heat Polished drainage eyes**  
Requires additional water soluble lubrication  
(For single use only)



**OptiSmooth<sup>®</sup> Mitro<sup>®</sup> XL<sup>®</sup>**  
Optimum smoothness Optimum comfort

Nelaton Catheter with 4 heat polished xtra large (XL) drainage eyes  
for Single Use Only

CAUTION: For use in USA see markings on device to sell by or on the order of a physician.



ISO 13485

Zhejiang Sunwood Technology Co., Ltd  
100% ISO 13485:2015 Certified  
100% ISO 13485:2015 Certified  
100% ISO 13485:2015 Certified

Corvus Ltd  
100% ISO 13485:2015 Certified  
100% ISO 13485:2015 Certified  
100% ISO 13485:2015 Certified

Vesica Urology Ltd

PO Box 4145, Slough,  
Berkshire, SL1 0QA  
[www.vesicaurology.com](http://www.vesicaurology.com)



***The New Mucous buster***

Four eyes are better than two  
(providing improved flow and drainage)




**Vesica Urology Ltd**

## OptiSmooth Mitro XL Key Features

- For the first time 4 unique oversized Xtra Large (XL) Heat Polished eyes. Developed especially for mitrofanoff bladder drainage. The oversized eyes will help to flush out the thick mucus produced when bowel tissue is used in the bladder reconstruction. Leading to potential improvements in mitrofanoff drainage and reduced infections.
- The four Xtra Large (XL) eyes are also heat polished and extremely smooth. Maximizing comfort on insertion and removal.
- Smooth uncoated PVC nelaton catheter (requires additional water soluble lubrication).
- Available in 40cm, 10ch, 12ch, 14ch (Charriere sizes)
- Flexible funnel connects securely to an external collection device such as a night drainage bag.
- Meaningful product codes OSMXL for OptiSmooth Mitro XL the next two numbers represent the catheter length and the final two numbers are the charrière size (catheter thickness).
- Can also be used for Ace procedures.

## OptiSmooth Mitro XL Product Range

40cm length

Size	Funnel Colour	Product Code	PIP Code
CH10		OSMXL4010	402-3768
CH12		OSMXL4012	402-3776
CH14		OSMXL4014	402-3784

Four perfectly formed extra large heat polished drainage eyes.



### OptiSmooth<sup>®</sup> Mitro<sup>®</sup> XL<sup>®</sup>

Optimum smoothness Optimum comfort

Nelaton Catheter with 4 heat polished extra large (XL) drainage eyes for Single Use Only  
CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

Zhejiang Songgood Technology Co., Ltd  
No. 200, Zhongyuan Road, Songgood Industrial Park, Songgood  
Zhejiang, China  
Phone: +86 571 8511 8888  
Fax: +86 571 8511 8889  
E-mail: sales@songgood.com

Vesica Urology Ltd  
PO Box 4145, Slough,  
Berkshire, SL1 0QA  
www.vesicaurology.com





**Vesica Urology Ltd**  
**PO Box 4145**  
**Slough**  
**SL10QA**

**Phone : +44 (0) 1753 445087**  
**Fax : +44 (0) 1753 447446**  
**Mobile : +44 (0) 7948 172607**  
**adrian@vesicaurology.com**  
**www.vesicaurology.com**

LATEX FREE

## Vesica Urology Ltd

OptiSmooth<sup>®</sup>, Mitro<sup>®</sup>, MitroXL<sup>®</sup> are registered trademarks  
 OPTISMOOTH MITRO IS A UK UNREGISTERED DESIGN 2011  
 OPTISMOOTH MITRO XL IS A UK UNREGISTERED DESIGN 2016